



| | |
|-----------------------|---|
| Product Description : | Face mask PP, blue, 3 plies, with ear loop |
| Product Code : | 15461932 (17-100-897) |
| Raw Material : | 28g/m ² + 20g/m ² Polypropylene |
| Size : | 17.5 x 9.5 cm |
| Packaging : | 50 pcs/bag, 10 bags/box |
| Carton Dimension : | 37 x 22 x 20 cm |



| Item | Length | Allowance | Remark |
|------|--------|-----------|--------|
| A | 17.5cm | +/- 1cm | |
| B | 9.5cm | +/- 1cm | |
| C | 11.5cm | +/- 1cm | |
| D | 16cm | +/- 1cm | |

© 2017 Thermo Fisher Scientific Inc. All rights reserved.
 Trademarks used are owned as indicated at fishersci.com/trademarks.

Austria: +43(0)800-20 88 40 Belgium: +32 (0)56 260 260 Denmark: +45 70 27 99 20
 Germany: +49 (0)2304 9325 Ireland: +353 (0)1 885 5854 Italy: +39 02 950 59 478
 Finland: +358 (0)9 8027 6280 France: +33 (0)3 88 67 14 14 Netherlands: +31 (0)20 487 70 00
 Norway: +47 22 95 59 59 Portugal: +351 21 425 33 50 Spain: +34 902 239 303
 Sweden: +46 31 352 32 00 Switzerland: +41 (0)56 618 41 11 UK: +44 (0)1509 555 500