

Burner, universal



For coal gas, natural (methane) gas, butane, propane or mixtures of these gases.

With adjustable jet, screwed air regulator, gas control needle valve, adjustable baffle. Manufactured in brass with nickel plating.



Catalogue No	Legacy No	O.D., mm	Height, mm
12852144	-	13 (burner tube), 9 (inlet tube)	160

Burner, bunsen



For natural (methane) gas.

Air regulator and rifled gas inlet tube. Nickel plated.



Catalogue No	Legacy No	O.D., mm	Height, mm	Pack qty
12832154	-	13 (burner tube), 9 (inlet tube)	145	5

Burner, Meker, brass



Designed for use with butane/propane mix gas cartridges.

This burner can achieve temperatures in excess of 1,000°C. The design of the gas cartridge enables the burner to be detached and re-attached as often as needed making it easy to store and transport. It is entirely constructed in brass with a replaceable plated steel perforated head. When used with the 220g cartridge the burner fits the standard 200mm (8 inch) high tripod. Typical burning time of the 220g cartridge is 2.75hr.



Catalogue No	Legacy No	Description
12802154	-	Burner complete with 220g cartridge
12812154	-	Gas cartridge, pack of 12

Burners, FIREBOY eco and FIREBOY plus

INTEGRA

Compatible with common gas types such as butane/propane and natural gas; can be connected to gas distribution system or gas cartridges.

FIREBOY eco

- Ignition by a footswitch (included) or push button
- Two adjustment knobs for air and gas allow easy fine-tuning of flame



FIREBOY plus

- Ignition by push button, or hands free operation with movement sensor or optional footswitch (not included)
- User defined maximal burning times
- Optional rechargeable battery
- Portable



The smooth, chrome plated metal housing is easy to clean and both UV and solvent resistant. With flame monitor and overheating protection. Functional and robust design. Convenient handling. Safe operation.

Catalogue No	Legacy No	Alt. No	Description
11833560	-	144000	FIREBOY plus
11843560	-	144010	FIREBOY eco

Accessories

Catalogue No	Legacy No	Alt. No	Description
11682080	-	144050	Adapter CV360
11813021	-	144051	Adapter CP250
11602090	-	144051	Adapter CV270/470
11863151	-	171081	Benchtop switch
10297601	-	143210	Borosilicate glass splash protector
11622100	-	144310	FIREBOY plus battery NiMH 3.6V

Burner, bypass



For natural (methane) gas.

As bunsen burner, with gas control and pilot ignition flame. Supplied with air regulator. Nickel plated.



Catalogue No	Legacy No	O.D., mm	Height, mm
12842154	-	13 (burner tube), 9 (inlet tube)	171