

## Tubes, centrifuge, Thermo Scientific Nalgene® Oak Ridge



Translucent polypropylene copolymer with polypropylene screw closure, which can be autoclaved repeatedly. Can be used up to 50,000xg in refrigerated and non-refrigerated centrifuges. Before autoclaving, just set the closure on top of the tube without engaging the thread. The actual brim capacity of the 50mL tube is 42mL.




Catalogue No	Legacy No	Alt. No	Capacity, mL	Pack qty
10706841	-	3119-0010	10	10
10613981	-	3119-0028	28	10
10757221	-	3119-0030	30	10
10757411	-	3119-0050	50	10

## Tubes, centrifuge, 15mL



- Black printed graduations and a large white marking area
- Certified non-pyrogenic and RNase/DNase free
- Sterile by gamma irradiation
- Available in racks or bulk packed
- Choice of cap styles - original plug seal or flat top cap



PET = Polyethylene terephthalate  
PP = Polypropylene

Catalogue No	Legacy No	Alt. No	Material	Cap style	RCF, x g [max.]	Pack type	Inner pack qty	Pack qty
10262861	-	430055	PET	Plug seal	3,600	Rack	50	500
10272861	-	430053	PET	Plug seal	3,600	Bag	25	500
10579691	-	430790	PP	CentriStar™ cap	12,000	Rack	50	500
10738771	-	430791	PP	CentriStar™ cap	12,000	Bag	25	500
10431013	-	430052	PP	Plug seal	12,000	Rack	50	500
10151151	-	430766	PP	Plug seal	12,000	Bag	25	500


### Accessories

Catalogue No	Legacy No	Alt. No	Description	Pack qty
10393023	-	431355	Standard rack, 15mL, polystyrene foam	20
10546813	-	4366	Universal rack, 15mL/50mL, polystyrene foam	20

## Tubes, centrifuge, 50mL



- Black printed graduations and large white marking spot
- Choice of cap styles
- Original plug seal or flat top rim seal
- Plug seal caps can be colour-coded using inserts
- Available in racks or bulk packed
- Certified non-pyrogenic and RNase/DNase free
- Sterile by gamma irradiation



PET = Polyethylene terephthalate  
PP = Polypropylene

Catalogue No	Legacy No	Alt. No	Material	Cap style	RCF, x g [max.]	Pack type	Inner pack qty	Pack qty
10334131	-	430828	PP	CentriStar™ cap	15,500	Rack	25	500
10604551	-	430829	PP	CentriStar™ cap	15,500	Bag	25	500
10314131	-	430290	PP	Plug seal	15,500	Rack	25	500
10694541	-	430291	PP	Plug seal	15,500	Bag	25	500
10684541	-	430304	PET	Plug seal	3,600	Rack	25	500
10291433	-	4558	PP	CentriStar™ cap	15,500	Universal rack	25	300

### Accessories

Catalogue No	Legacy No	Alt. No	Description	Pack qty
10164212	-	4365	Standard rack, 50mL, polystyrene foam	20
10546813	-	4366	Universal rack, 15mL/50mL, polystyrene foam	20