

# Centrifuge tubes



For additional labels see pages 914 to 915

## Labels, microcentrifuge tube, Tough-Tags® and Tough-Spots®



Labels for microcentrifuge tubes and other containers in the research and clinical laboratory.

Unlike conventional tape materials, these labels withstand steam autoclaving, boiling water baths, freezer temperatures, organic solvents, caustic agents and other challenges without peeling.

Pre-cut, peel off labels perfectly sized to fit microcentrifuge tubes and other containers in the research and clinical laboratory. Chemically inert polyester with strong adhesion to plastic and other materials. Microcentrifuge tubes with Tough-Tags® wrap around labels slide easily in and out of rotors without binding because labels are thin and the adhesive does not migrate even after boiling and autoclaving. Available in laser sheets for labelling of large columns of samples direct from a laser printer.

### Tough-Tags®, rolls

Catalogue No	Legacy No	Description	Colour	Dimensions, mm	Pack qty
12338807	-	Small	Clear	25 x 12.7	1,000
12348807	-	Small	White	25 x 12.7	1,000
12378807	-	Small	Red	25 x 12.7	1,000
12388807	-	Small	Yellow	25 x 12.7	1,000
12398807	-	Small	Green	25 x 12.7	1,000
12308817	-	Small	Blue	25 x 12.7	1,000
12358807	-	Large	Clear	32.5 x 12.7	1,000
12368807	-	Large	White	32.5 x 12.7	1,000
12318817	-	Large	Red	32.5 x 12.7	1,000
12328817	-	Large	Yellow	32.5 x 12.7	1,000
12338817	-	Large	Green	32.5 x 12.7	1,000
12348817	-	Large	Blue	32.5 x 12.7	1,000

### Tough-Spots®, rolls

Catalogue No	Legacy No	Description	Colour	Dimensions, mm	Pack qty
12358817	-	Small	White	9.5	1,000
12368817	-	Small	Red	9.5	1,000
12378817	-	Small	Yellow	9.5	1,000
12388817	-	Small	Green	9.5	1,000
12398817	-	Small	Blue	9.5	1,000
12308827	-	Large	White	12.7	1,000
12318827	-	Large	Red	12.7	1,000
12328827	-	Large	Yellow	12.7	1,000
12338827	-	Large	Green	12.7	1,000
12348827	-	Large	Blue	12.7	1,000

### Cryogenic labels for laser printers, sheets

Catalogue No	Legacy No	Description	Colour	Dimensions, mm	Pack qty
12672580	-	Tough Spots™ large	White	9.5	3,840
12662580	-	Tough-Tags™ small	White	24 x 12.7	2,975
12682580	-	Tough-Tags™ large	White	32.5 x 12.7	2,125



## Label station, Tough-Tags®, Fisherbrand®

**Fisherbrand®**

Portable Tough-Tags® station holds up to six standard size boxes of Tough-Tags®, Cryo-Babies®, Tough-Spots® and TeenyTough-Spots®. Just drop the box into the station and feed the paper through the slot. The acrylic station includes a writing surface, cutting edge and non-slip feet for easy dispensing. Counter sunk holes allow the unit to be wall mounted if required.



Tough-Tags®, Cryo-Babies®, Tough-Spots® and TeenyTough-Spots® are registered trademarks of Diversified Biotech. Tags sold separately to holder and shown for display purposes.

Catalogue No	Legacy No	Description
11947994	-	Tough-Tags® station