

Tube racks, centrifuge, Rota-Rack™, Fisherbrand®



Fisherbrand®

Each module can be independently rotated enabling the use of several sizes of test tubes, culture tubes, centrifuge tubes or microcentrifuge tubes at the same time. Click each module in to place, load with tubes and carry to a new location. The racks remain fully orientated until they are twisted again.



Each module of the large Rota-Rack™ holds 2 x 50mL tubes, 8 x 15mL tubes, 6 x 16mm to 20mm tubes, or 10 x 10mm to 13mm tubes.

Each module of the small Rota-Rack™ holds 6 x 15mL tubes, 9 x 1.5/2mL tubes, 12 x 0.5/0.6mL tubes or 32 independent 0.2mL PCR* tubes or 4 x 8 tube strips.

Each rack has modules in green, pink, blue and yellow, and is fully autoclavable.

Test tube racks, Rota-Rack®

Catalogue No	Legacy No	Description	Dimensions, internal [w x d x h], mm
11384085	-	Large	224 x 110 x 110
11394085	-	Small	210 x 100 x 100

Tube racks, centrifuge, Puzzle Rack, Fisherbrand®



Fisherbrand®

- For use with 1.5mL, 2.0mL, 15mL and 50mL microtubes and centrifuge tubes
- Use individually, configure horizontally or in the shape of a square
- Racks easily connect to each other and can expand with additional puzzle racks
- Tubes have references printed on to aid retrieval
- Bright and colourful design
- Racks can be stacked when not in use
- Polypropylene
- Dimensions: 151mm x 135mm x 66mm



Two-sided puzzle shaped pieces hold an assortment of test tubes.


Catalogue No	Legacy No	Description	Pack qty
11324095	-	Puzzle Rack	4

Tube rack, PCR*, Rota-Rack Duo, Fisherbrand®



Fisherbrand®

Rotating rack holds 48 x 0.2mL, 0.5mL and 1.5mL PCR tubes. Polypropylene with imprinted alphanumeric grid aid. 296mm x 118mm x 34mm.



Catalogue No	Legacy No	Description
11304095	-	Rota-Rack Duo rotating tube rack

*Polymerase Chain Reaction (PCR) is a process covered by patents owned by Hoffmann-La Roche

Tube racks, microcentrifuge, polypropylene, 100 well, Fisherbrand®



Fisherbrand®

Durable 100 well storage racks feature three point hinges and clasp for secure lid closure.

Hold 1.5mL to 2mL microtubes and cryovials. User friendly racks feature moulded grid lines on the lid, imprinted coordinates on the bottom of the base and on the box. Racks measure 141mm x 151mm x 57mm. Autoclavable.



Catalogue No	Legacy No	Description	Pack qty
11790334	-	Natural	5
11918014	-	Black	5
11700344	-	Assorted colours (blue, green, purple, yellow and orange)	5