

# Containers storage



## Containers, storage box, polypropylene, 0.11L to 13L

A range of microwaveable, freezer safe, stain resistant storage boxes with snap on lids. All manufactured in the UK.

### Round

Catalogue No	Legacy No	Dimensions [w x d], mm	Capacity, mL
12378288	-	80 x 35	110
12388288	-	80 x 75	250
11368053	-	80 x 95	370
11378053	-	230 x 95	3,000
11388053	-	205 x 150	4,350

### Square

Catalogue No	Legacy No	Dimensions [w x d x h], mm	Capacity, mL
11398053	-	130 x 130 x 50	500
11308063	-	135 x 135 x 140	1,800
10129961	-	250 x 250 x 75	3,500
11328063	-	260 x 260 x 125	6,500
11338063	-	315 x 315 x 160	13,000

### Rectangle

Catalogue No	Legacy No	Dimensions [w x d x h], mm	Capacity, mL
11348063	-	175 x 120 x 60	750
11358063	-	195 x 140 x 60	1,000
11368063	-	200 x 120 x 70	1,100
11378063	-	235 x 170 x 80	2,250
11388063	-	270 x 195 x 105	3,750
11398063	-	300 x 215 x 140	7,500

### Beaker with screw lid

Catalogue No	Legacy No	Dimensions [w x d], mm	Capacity, mL
11308073	-	78 x 145	330



## Containers, storage box with lid, polypropylene, 15L to 100L

Transparent storage container with lid, manufactured from polypropylene. Stackable.

Catalogue No	Legacy No	Alt. No	Dimensions [w x d x h], mm	Capacity, L
11346644	-	1900/M/TRAS/15	350 x 250 x 200	15
12861834	-	1900/M/TRAS/20	400 x 300 x 220	20
12871834	-	1900/M/TRAS/30	450 x 350 x 250	30
12881834	-	1900/M/TRAS/40	500 x 400 x 260	40
12891834	-	1900/M/TRAS/50	550 x 400 x 290	50
12801844	-	1900/M/TRAS/60	620 x 400 x 310	60
12811844	-	1900/M/TRAS/100	750 x 440 x 400	100



## Containers, transport, polyethylene, 50kg to 100kg, Strong/Duraflex

Strong or Duraflex transport containers, with polypropylene lids available separately.

### Strong container

Catalogue No	Legacy No	Alt. No	Dimensions [w x d x h], mm	Capacity, kg
12821844	-	8325/HD/50	500 x 400 x 300	50
12831844	-	8325/HD/70	600 x 500 x 360	70
12841844	-	8325/HD/100	750 x 500 x 420	100

### Duraflex container

Catalogue No	Legacy No	Alt. No	Dimensions [w x d x h], mm	Capacity, kg
12871844	-	8350/HD/100	750 x 500 x 420	100

### Lids

Catalogue No	Legacy No	Alt. No	Dimensions [w x d], mm
12881844	-	9175/M/50	500 x 400
12891844	-	9175/M/70	600 x 500
12801854	-	9175/M/100	750 x 500

