

## Containers, storage case, polyethylene, stackable, 10L to 60L

Stackable storage cases manufactured from polyethylene. Polypropylene lids are available separately.

### Containers

Catalogue No	Legacy No	Alt. No	Dimensions [w x d x h], mm	Capacity, L
12821854	-	1885/HD/10	300 x 400 x 100	10
12831854	-	1885/HD/12	600 x 400 x 70	12
12841854	-	1885/HD/20	600 x 400 x 100	20
12454060	-	1885/HD/28	600 x 400 x 130	28
12861854	-	1885/HD/45	600 x 400 x 190	45
12871854	-	1885/HD/50	600 x 400 x 240	50
12881854	-	1885/HD/60	600 x 400 x 320	60

### Lids

Catalogue No	Legacy No	Alt. No	Dimensions, mm
12891854	-	9195/M/3040	300 x 400
12388268	-	9195/M/6040	600 x 400



## Containers, storage, Fisherbrand® Tubby®

**Fisherbrand®**

- Stackable storage containers to help with laboratory organisation or transportation of products
- Removable dividers
- For keeping gloves, tubes, tips, pipettes and other small items organised and tidy
- Stackable and proportionately sized to utilise space efficiently



Catalogue No	Legacy No	Description	Dimensions [w x d x h], mm	Pack qty
11938014	-	Tubby®, with 5 tubs, 5 lids and 15 dividers	330 x 200 x 115	5
11948014	-	Tubby® Jr., shallow, with 5 tubs, 5 lids and 15 dividers	330 x 200 x 70	5
11918004	-	Tubby® Sr., large, with 2 tubs, 2 lids and 12 dividers	300 x 415 x 115	2

## Trays, utility, Fisherbrand®



**Fisherbrand®**

- Lightweight solution for transporting tubes and smaller items efficiently

The cover for utility tray protects its contents and provides access to the tray handle in its vertical position. Two holes on the lid top corners allow cable ties to thread through for secure storage. Measures 420mm x 240mm x 205mm. Extra pockets allow for extra storage on the tray. They simply clip on to both ends of the tray. Measures 230mm x 110mm x 105mm. Urine/blood collection inserts hold up to 11 sample containers, 9 inserts for 59mm dia. containers and 2 square spaces for 64mm dia. containers. Insert measures 368 x 210 x 80mm.

Catalogue No	Legacy No	Description
11740634	-	Utility tray with 13mm OneRack™
11750634	-	Utility tray with 16mm OneRack™

### Accessories

Catalogue No	Legacy No	Description	Pack qty
11760634	-	Disposable inserts for 13mm OneRack™	4
11770634	-	Disposable inserts for 16mm OneRack™	4
11780634	-	Lid for blood collection tray, polystyrene, stackable, 420mm x 240mm x 205mm	1
11790634	-	Extra storage for blood collection tray, acrylonitrile butadiene styrene, 230mm x 110mm x 105mm	2
11700644	-	Insert for blood collection tray, holds 11 urine samples, 9 round and 2 square, comes with grid inserts for square spaces, polystyrene, 368mm x 210mm x 80mm	4

