

Refrigerators, universal, spark free



- Static cooling system
- ATEX certified
- Seam free plastic interior, so easy to clean and keep hygienic
- Sturdy lock on door protects the contents from unauthorised access
- Replaceable door seal
- Adjustable feet
- Temperature range +2°C to +10°C
- Environmentally friendly hydrocarbon refrigerant

Technical Specification - Specific

	12978604	12028359	12078261
Capacity, L	260	360	500
Shelves supplied	4	5	5
Dimensions, external [w x d x h], mm	600 x 600 x 1,216	600 x 600 x 1,590	752 x 710 x 1,516
Dimensions [w x d x h], mm	513 x 441 x 1,047	513 x 441 x 1,418	652 x 537 x 1,338
Power rating, W	140	150	150
Energy consumption, kWh/24hr	1.00	1.10	1.10

Catalogue No	Legacy No	Alt. No	Capacity, L
12978604	-	FKEX 2600	260
12028359	-	FKEX 3600	360
12078261	-	FKEX 5000	500



Technical Specification - General

Temperature range, °C	+2 to +10
Electrical supply	230V 50/60Hz

Refrigerators, laboratory, upright, fan assisted



- Large capacity
- Legs as standard
- Fitted lock as standard
- Environmentally friendly hydrocarbon refrigerant

A choice of white or stainless steel exterior refrigerators, available in two sizes.

Technical Specification - Specific

	12078281	11916239	12088281	11926239
Capacity, L	583	583	663	663
Dimensions, external [w x d x h], mm	750 x 750 x 1,864	750 x 750 x 1,864	750 x 750 x 1,864	750 x 750 x 1,864
Dimensions, internal [w x d x h], mm	650 x 581 x 1,510	650 x 581 x 1,510	650 x 581 x 1,710	650 x 581 x 1,710
Energy consumption, kWh/24hr	1,269	1,269	1,298	1,298
No. of shelves	5	5	5 plus base grid	5 plus base grid
Exterior finish	White	Stainless steel	White	Stainless steel
Door finish	White	Stainless steel	White	Stainless steel
Mass, kg	88	87	97	95

Technical Specification - General

Cooling system type	Fan assisted
Defrosting system	Automatic
Temperature range, °C	+1 to +15
Interior finish	Moulded
Temperature display type	External digital
Castors	Standard
Door hinges	Right, reversible on site
Electrical supply	220-240V, 50Hz

Catalogue No	Legacy No	Alt. No	Description	Capacity, L
12078281	-	GKv5730	Refrigerator, fan assisted, white exterior	583
11916239	-	GKv5790	Refrigerator, fan assisted, stainless steel exterior	583
12088281	-	GKv6410	Refrigerator, fan assisted, white exterior	663
11926239	-	GKv6460	Refrigerator, fan assisted, stainless steel exterior	663

