

Refrigerators, laboratory, fan assisted



- Data logging of temperature every two minutes
- Audible and visual door-open and hi-low temperature alarm
- Volt free contact for connection to remote system
- Data transferable to external unit

A range of purpose built laboratory refrigerators, single and double door, solid and glass door.



Technical Specification - Specific

	10646875	10207694
Capacity, L	601	601
Dimensions, external (w x d x h), mm	700 x 830 x 2,150	700 x 830 x 2,150
Energy consumption, kWh/24hr	2.4	2.8
Temperature range, °C	-2 to +16	0 to +16
No. of shelves	4	4
Mass, kg	155 (gross)	155 (gross)
Door finish	Solid	Glass

	10094084	10145924
Capacity, L	1,427	1,427
Dimensions, external (w x d x h), mm	1,430 x 830 x 2,150	1,430 x 830 x 2,150
Energy consumption, kWh/24hr	3.4	4.1
Temperature range, °C	-2 to +16	0 to +16
No. of shelves	8	8
Mass, kg	259 (gross)	250 (gross)
Door finish	Solid	Glass

Technical Specification - General

Defrosting system	Automatic
Exterior finish	White steel
Interior finish	304 stainless steel
Castors	Standard
Temperature display type	Digital
Door lock	Standard
Door hinges	Right, reversible on site

For Fisherbrand® unisex laboratory coats, please see page 1305.



Catalogue No	Legacy No	Alt. No	Description
10646875	-	LKPv6520	Refrigerator, single solid door
10207694	-	LKPv6522	Refrigerator, single glass door
10094084	-	LKPv1420	Refrigerator, double solid door
10145924	-	LKPv1422	Refrigerator, double glass door

FINDING THE PRODUCTS YOU NEED