


Tubing, peristaltic pump, PVC, for 205U/BA cassette pump series



Catalogue No	Legacy No	Bore size, mm	Colour	Length, mm	Pack qty
12406310	-	0.13	Orange/black	400	12
12416310	-	0.25	Orange/blue	400	12
12426310	-	0.50	Orange/yellow	400	12
12436310	-	0.76	Black/black	400	12
12446310	-	1.02	White/white	400	12
12456310	-	1.29	Grey/grey	400	12
12466310	-	1.42	Yellow/yellow	400	12
12476310	-	1.65	Blue/blue	400	12
12486310	-	2.05	Purple/purple	400	12
12496310	-	2.54	Purple/orange	400	12
12406320	-	2.79	Purple/white	400	12



Tubing, peristaltic pump, Marprene



Marprene is a long life tubing, highly resistant to attack by oxidising agents including ozone, peroxides and hypochlorite. It has high acid/alkali resistance and good resistance to animal and vegetable oils. Autoclavable, Marprene may be steam sterilised at 135°C or cold sterilised by ethylene oxide.

Tubing for 323, 400 and 500 series pumps

Catalogue No	Legacy No	Bore size, mm	Wall thickness, mm	Length, m
10627762	-	0.5	1.6	15
12426270	-	0.8	1.6	15
12436270	-	1.6	1.6	15
10488501	-	3.2	1.6	15
10027911	-	4.8	1.6	15
12466270	-	6.4	1.6	15
12476270	-	8.0	1.6	15

Tubing for 505Di/L pump

Catalogue No	Legacy No	Bore size, mm	Wall thickness, mm	Length, m
12436330	-	1.6	2.4	15
12446330	-	3.2	2.4	15
12456330	-	4.8	2.4	15
12466330	-	6.4	2.4	15
12476330	-	8.0	2.4	15
12486330	-	9.6	2.4	15

Marprene double 'y' elements

Catalogue No	Legacy No	Bore size, mm	Wall thickness, mm	Pack qty
12496330	-	1.6	2.4	5
12406340	-	3.2	2.4	5
12416340	-	4.8	2.4	5
12426340	-	6.4	2.4	5
12436340	-	8.0	2.4	5
12446340	-	9.6	2.4	5

Silicone double 'y' elements

Catalogue No	Legacy No	Bore size, mm	Wall thickness, mm	Pack qty
12486320	-	1.6	2.4	5
12496320	-	3.2	2.4	5
12406330	-	4.8	2.4	5
12416330	-	8.0	2.4	5
12426330	-	9.6	2.4	5

