



Aluminium Fold-Up Cart


Sturdy aluminium alloy construction permits this lightweight cart to be used for loads up to 150kg pounds. Move your heaviest lab chemicals, drums, and cases of glassware with this cart.

Adjustable height

Folds up for easy storage to a compact briefcase size. Handle adjusts to two heights (68.6 cm to 83.8 cm) to accommodate all users. Telescoping platform extends for longer loads from 45.7 to 64.8 cm. Width is 406 mm.

Supplied bungee cords keep items from slipping. Large, solid rubber wheels make this cart extremely easy to push. Rubber bumpers on all four corners protect walls from accidental nicks. Cart weight is just 8.2kg.

Description	Cat. No.	Price
Fisher Scientific Aluminium Fold-Up Truck	FB61149	£145.00



Aluminium Fold-Up Hand Truck

Move cartons, chemicals, and gas cylinders easily. Aluminium alloy construction allows this lightweight unit (only 6.8kg) to move loads up to 100kg. Supplied bungee cords hold materials firmly in place.

Stores in a drawer

Folds to a compact size for easy storage. Height of handle is 71.1 cm (lower position) and 106.7 cm (extended). Flip-down platform size is 35.6 cm wide by 24.1 cm. Big, smooth-rolling, solid-rubber 16.5 cm diameter tires make it easy to move heavy loads.

Description	Cat. No.	Price
Fisher Scientific Aluminium Fold-Up Hand Truck	FB61150	£128.00



Mini-Dolly Folding Cart

This is the ideal platform cart to move chemicals, carboys, computer paper, filing boxes, and scientific equipment. Eliminates lifting required by high-shelf carts. Cart ingeniously folds to a compact size for easy storage.

Easy rolling

Weighs only 9kg, but handles up to 120kg. Large nonskid platform (45.1 cm x 64.8 cm) holds materials firmly. Plastic guard on four sides protects your walls. Supplied elastic cords with quick hooks secure the load.

Super-glide wheels swivel for complete manoeuvrability. Easy-grip handle is 87.6 cm long and ensures a firm hold. Constructed of high-impact ABS plastic and chrome-plated steel for durability.

Description	Cat. No.	Price
Fisher Scientific Mini-Dolly Folding Cart	FB61148	£115.00