



### Storage racks, Cryogenic, microplate, aluminium, Thermo Scientific Nunc®



76

For storage in chest and upright mechanical freezers or incubators.

- Adjustable shelves to fit any plate or tube rack
- Colour coded top and sides for easy identification
- Locking bar keeps plates secure
- Autoclavable
- Accommodates: 18 x low profile, 10 x standard, 6 x DeepWell™ 1mL, 5 x DeepWell™ 2mL
- Dimensions (w x d x h), mm: 9.5 x 13.5 x 23.5

Catalogue No	Alt. No	Colour
<b>MPA-922-010Q</b>	367001	Blue
<b>MPA-922-020N</b>	367002	Natural
<b>MPA-922-030K</b>	367003	Yellow
<b>MPA-922-040H</b>	367004	Red



### Storage racks, vial, cryogenic, aluminium, Thermo Scientific Nunc® CryoRacks



76

For storage of CryoBoxes™ and CryoTubes™

- Adjustable shelves to fit many box sizes
- Accommodates up to 12 Max -100 boxes for 1.8mL tube or 8 MegaMax -100 boxes for 3.6mL tubes, or a combination of these
- Alignment rods secure boxes in position
- Easy to clean
- Fit most chest freezers
- 11 adjustable shelves and 12 compartments
- Lightweight
- Dimensions (w x d x h), mm: 143.2 x 140.4 x 670

Catalogue No	Alt. No	Plate colour
<b>MPA-921-010W</b>	367000	Natural (top)



### Storage racks, vial, cryogenic, Thermo Scientific Nalgene®, CryoBox™



55

- Stainless steel racks, which hold CryoBoxes™ securely but separately
- Individual boxes can be removed and the rack returned to the freezer without disturbing other boxes
- Horizontal racks fit most upright freezer compartments

Catalogue No	Alt. No	No. of shelves	For box [array]	Dimensions [w x d x h], mm
<b>CRY-200-010C</b>	DS5035-0004	4	5 x 5	84 x 83 x 224
<b>CRY-200-030T</b>	5036-0004	4	9 x 9	143 x 140 x 224
<b>CRY-200-040Q</b>	5036-0009	9	9 x 9	143 x 140 x 502
<b>CRY-210-011Q</b>	5038-4322	12 (3 x 4)	9 x 9	445 x 143 x 225
<b>CRY-210-021N</b>	5038-4422	16 (4 x 4)	9 x 9	590 x 143 x 225

